



# TeatreeVag

Intimate gel with tea tree oil

20ml

**Ingredients (INCI):** Aqua, Poloxamer 407, Propylene Glycol, Melaleuca Alternifolia, Panthenol, Lactic Acid, Sodium Benzoate, Thymol, Limonene.

**Action:** **Teatreevag** intimate gel contains Tea tree oil, thymol and panthenol in bio adhesive gel - base. Tea tree oil (lat. *Melaleuca alternifolia*) performs significant anti-inflammatory, antimicrobial and epithelializing effect. Tea tree oil's immune stimulative effect provides natural vaginal flora balance, and natural vaginal acidity, that protects it from occurrence of bacterial, viral and fungal infections. Thymol provides antimicrobial properties, and it helps resorption of tea tree active ingredients into the inflamed tissue. Panthenol regenerates irritated and damaged mucosa and it keeps its natural humidity. Bio adhesive gel excipients have an epithelizing effect, and it provides continuous, step by step gel's active ingredients release. Applied on intimate region, gel coats the mucous membrane, regenerates it and eliminates all the inconveniences that occur as a result of inflammation and infection: itching, burning, dryness, redness, Odor.




**Teatreevag** intimate gel is recommended for:

- as an aid in the inflammatory process of the vaginal mucous membranes
- as an aid in the treatment of infections of the vaginal mucosa accompanied by vaginal secretion
- in all conditions when the vaginal balance is disturbed
- prevention in persons prone to infection intimate regions

**Application:** Apply vaginally, in the evening before going to bed, in a lying position. After opening the tube and winding the applicator, squeeze 2 cm of gel (dose for one application) deep into the vagina with gentle pressure. The amount of gel is enough for 5 applications. With superficial wounds apply gel in thin layer 2-3 times daily.

**Notes:** Do not use in case if you are sensitive (allergic) on any of ingredient. Not recommended for pregnant women. For external use only.

**Package:** Box with 20 ml gel tube and two applicators.



**Storage instructions:**

Store at temperatures up to 25°C, out of reach of children.



**Shelf life:** 30 months.



CORNUS GmbH  
Zentgrafenstrasse 128,  
34130 Kassel, Deutschland

**CORNUS**<sup>®</sup>  
natural nutrition