

Stop biting nails GEL



Gel against nail biting and for nail care, 5g For both child and adult nails

Active ingredients:	100g of product contain:	Content (mg) per 0.15 g product
Aqua purificata	52.202 g	0,783 mg
Carbomer	0,5 g	0,00075 mg
Glycerin	35 g	0,053 mg
Gentiana Lutea Root Extract, Aqua, Glycerin	1 g	0,0015 mg
Panthenol	3 g	0,0045 mg
Propoils Extract	5 g	0,0075 mg
Propylene Glycol	2 g	0,003 mg

Ingredients: Aqua, Gentiana Lutea Extract, Poloxamer 407, Propolis Extract, Panthenol, Propylparaben, Methylparaben.

Effects:

Stop Biting Nails Gel is a gel that deters biting of the nails, cuticles as well as thumb sucking.

It is applied to the surface of nails and the skin surrounding nails. Due to its bitter taste stemming from the extract of the Gentiana lutea, **Stop Biting Nails Gel** eliminates the bad habit of nail biting, which is not only an aesthetic, but also a health issue (nail biting can cause a fungal infection of the fingers and in the mouth).

Panthenol contributes to the regeneration of sensitive skin surrounding the nails.

Stop Biting Nails Gel contains a specially formulated gel as base, which enhances the contact of the product with skin and nails and gradually releases active ingredients, thus enabling an extended effect of the gel on the application area. **Stop Biting Nails Gel** is a natural and completely safe

product, recommended both for children and adults.

The propolis extract in the product nourishes skin making it

soft and elastic.

Regular application of the gel can prevent the breaking and

peeling of nails, as well as formation of cuticles.

They are recommended as help:

- Against nail biting, cuticle sucking and thumb sucking.
- ✓ For beautiful and well-groomed nails and healthy growth.

Instruction for use:

A thin layer of gel should be applied onto the surface of the nail and the skin surrounding the nail several times during the day.

Warnings:

For external use only. Not recommended for use by those hypersensitive to any of the ingredients in the product.

Storage:

Store the closed tube up to 30°C.

Shelf life:

4 years from the date of manufacture.

The date of manufacture is contained in the serial number imprinted on the weld of the tube and printed on the box.

Packaging:

Tube containing 5g of gel



CORNUS GmbH, Zentgrafenstraße 128, 34130 Kassel, Deutschland